

## **HOSPITAL PATIENT IDENTIFICATION**

### **Situation**

A Regional Hospital in Central Illinois was using Laser wristbands provided from PDC. They chose PDC because they were inexpensive. They were having many issues with jams in laser printer and print quality, causing bar codes to be unreadable.

### **Solution**

DFI is a Platinum Distributor for LaserBand.

DFI's equipment specialist recommended converting to Direct Thermal 300 DPI TSC Printers and utilizing Laser Band's "Fusion Band" which is a self laminating Direct Thermal wristband. The lamination provides protection from direct exposure to lotions, baths etc that could compromise the quality of the wristband. Also utilizing Nicelabel/Nicewatch software, a command is sent to the appropriate printer to print a Pediatric or Adult wristband.

### **Results**

The hospital has eliminated all jams and provides a legible wristband that can be read in all cases.

### **Customer Comments**

"The clarity of the bands is unbelievable" stated the unit clerk